

RAISING THE BAR WITH PATIENT OUTCOMES

Full Range of Focus
Non-Diffractive IOLs




First IOLs mimicking the optical
performance of the human lens

VOPTICA
SMART VISUAL OPTICS



The bilateral implantation of ArtIOLs® 70 together with ArtIOLs® 40 delivers a full range of focus and excellent VA at all distances¹.



Monocular distance-corrected contrast sensitivity for ArtIOLs® 70 & ArtIOLs® 40 with different depth of focus, were identical at all spatial frequencies² and comparable to monofocal IOLs³


1. Marín JM, et al. Visual Performance at All Distances and Patient Satisfaction With a New Aspheric Inverted Meniscus Intraocular Lens. J Refract Surg. 2023 Sep;39(9):582-588.

2. Data onfile PRJ-INV_005 (37 patients)

3. Macias-Murelaga B, et al. Prediction Error Stabilization and Long-Term Standard Results with a Monofocal Intraocular Lens. Vision 2022;6(1):5.

HIGH SPECTACLE INDEPENDENCE

Percentage of patients reporting
spectacle independence²



Patient-Reported Spectacle Independence Questionnaire (PRSIQ)⁴

2. Data on file PRJ-INV_005 (43 patients)


4. Morlock R, et al. Patient-Reported Spectacle Independence Questionnaire (PRSIQ): Development and Validation. Am J Ophthalmol. 2017 Jun;178:101-114.

HIGH QUALITY OF VISION


Percentage of patients never
having experienced²



HALOES




GLARE




STARBURSTS

The Quality of Vision (QoV) Questionnaire⁵


ArtIOLs® non-interrupted
and non-stepped optical surface
minimizes visual disturbances.



Haloes



Glare




Starbursts

2. Data on file PRJ-INV_005 (43 patients)


5. McAlinden C, et al. The Development of an Instrument to Measure Quality of Vision: The Quality of Vision (QoV) Questionnaire. Invest. Ophthalmol. Vis. Sci. 2010;51(11):5537-5545

ADVANCED OPTICAL TECHNOLOGY


The ArtIOLs® reproduce the imaging formation properties of the **natural crystalline lens** in the periphery of the retina with their **inverted meniscus design**⁶




Standard IOL
(Biconvex)



ArtIOLs®
(Inverted meniscus)




Biconvex IOL



ArtIOLs®
Inverted Meniscus IOL

Patients implanted with ArtIOLs® present a **reduction of spherical equivalent by 1.97D at 40° of temporal eccentricity** compared to standard biconvex IOLs⁶

The ArtIOLs® inverted meniscus design minimizes negative dysphotopsias⁷



Standard IOL (Biconvex)⁸



ArtIOLs® (Inverted meniscus)

1. Marín JM, et al. Visual Performance at All Distances and Patient Satisfaction With a New Aspheric Inverted Meniscus Intraocular Lens. J Refract Surg. 2023 Sep;39(9):582-588.

6. Villegas EA et al. Peripheral, Refraction and Contrast Detection Sensitivity in Pseudophakic Patients Implanted With a New Meniscus Intraocular Lens. J Refract Surg. 2022 Apr;38(4):229-234

7. Harilaos G, et al. Quantification And Analysis Of Negative Dysphotopsia And Associated Phenomena In A Physical Model Of The Eye. Invest. Ophthalmol. Vis. Sci. 2022;63(7):1798 – F0414.

8. Pusnik A, et al. Dysphotopsias or Unwanted Visual Phenomena after Cataract Surgery. Life (Basel). 2022 Dec 24;13(1):53

ArtIOLs®

TECHNICAL SPECIFICATIONS



voptica.com

Support:

info@voptica.com

(0034) 868 881 714

Optics

Lens type	Single piece foldable lens
Optical design	Aspheric optics with extended depth of focus
Shape	Inverted meniscus
Material	Hydrophobic Acrylic UV absorbing and blue light filter
Power range	+12.00 to +30.00 Diopters (0.50D steps)
Optical diameter	6.0 mm
Total diameter	13.0 mm
Refractive index	1.54 (glistening free)
Edge design	Square

Optical Biometry

Suggested A-Constant*	ArtIOLs® 70 / ArtIOLs® 40
SRK/T	120.0
Haigis	$a_0 = 0.720, a_1 = 0.297$ and $a_2 = 0.169$

Haptics

Haptic design	C-L haptic
---------------	------------

Delivery System

Injector-Cartridge set	Single use
Recommended incision size	$\geq 2.2\text{mm}$
Lens delivery	Single handed plunger

*It is recommended that surgeons personalize the constants they use.



2460
Caregroup Sight
Solution PL

MP-TR-0010

VOPTICA
SMART VISUAL OPTICS